

Last Update: April 10, 2023

# **Rolling Circle Amplification ( RCA )**


# CONTENTS

1. RCA concept
2. Experiment in Bulk solution
3. Characterization of RCA product
4. RCA in GUV
5. Notes

# RCA Concept

## Introduction:

- RCA stands for Rolling Circle Amplification which is a process of generating copies on DNA that are complementary to the template.
- This is similar to a PCR reaction in which amplification of DNA is done using a template, staple in the presence of free bases and an enzyme which does the polymerization.
- There can be many cross polymerizations that might occur in PCR reaction. These might lead to undesired DNA strands being generated than desired.
- RCA employs a circular template which avoids most of the chance of cross polymerizations.



# Experiment in Bulk solution

## Initial Materials:


**SC-DT1**, 46nt: GCCTTCCCTAACCCTAACCCTAACCCTAACCCTAACCCTAACCCTAACCCTAAATGACT (Template Strand)

**SC-DP1**, 18nt: AGGGAAGGC AGTCATTAG (Splint Strand)

## Step 1: Phosphorylation of 5' end



- Can be done by using Kinase Enzyme.
  - Kinase enzyme should be used with the Ligase buffer.
  - Linear Strand is introduced such that it is 20  $\mu\text{M}$  in the total volume (50 $\mu\text{L}$ ).
  - Anneal it at 37°C for 2 hr.
- 
- T4 Kinase Enzyme: 10000 units/mL

Note: Using Ligase Buffer from first step instead of kinase buffer is recommended. Ligase buffer contains ATPs and DTT's. Kinetion and Ligation can happen in Ligase buffer but, Ligation cannot occur in Kinase buffer alone.





## Step 2: Binding Splint Strand with Circular Strand

- Introduce the Staple strand (1.2 times molar ratio of Linear Strand)
- Put it in a annealing cycle
  - 95°C for 5 min
  - 65°C for 30 min
  - 50°C for 30 min
  - 37°C for 30 min
  - 22°C for 20 min
- Staple strand will bind to the Phosphorylated Linear strand to form a circular shape.




## Step 3: Ligating The Circular Strand

- This step bonds the phosphorylated 5' end to the 3' end and makes the circular strand to form a complete circular template.
- Ligase Enzyme is used for this job.
- Add Ligase Enzyme in required amount and leave the sample in a dark place for 16hrs.
- Ligase: 400000 units/mL




## What could possibly be in the final solution?

- ✓ Free DNA (Linear and Splint Strands)
- ✓ Non-Ligated Circular Template
- ✓ Secondary Structures
- ✓ Ligated Circular Template ( Desired )



## Step 4: Extracting Ligated Circular Template

- Use dPAGE to extract the circular template.
- Follow the dPAGE protocol and electro dialysis for recovery.
- Butanol concentration is preferred to extract the DNA from dialysis product.
- Dissolve the DNA salt in DI water and measure the concentration by absorption measurements.



## RCA calculations:

- Look into the excel sheet for any change in new concentrations

### Stock Solutions:

	<i>Linear Strand ( SC-DT1 )</i>	<i>Staple Strand (SC-DT2 )</i>
<b>OD</b>	34.9	17
<b>nmoles</b>	83.9	89.6
<b>mg</b>	1.16	0.5
<b>Mol. Wt. ( g/mol )</b>	13783	5612.7
<b>Conc. (mM)</b>	1	0.1
<b>Conc. ( ug/uL)</b>	13.826	

<b>Kinase</b>	10000	units/mL
<b>Ligase</b>	400000	units/mL
<b>Phi29</b>	10000	units/mL
<b>Buffers</b>	10	X
<b>dNTPs</b>	10	mM

### Phosphorylation:

	<i>Volume (uL)</i>	<i>Conc. (uM)</i>
<b>Linear Strand</b>	1	20
<b>Kinase</b>	3	
<b>Ligase Buffer</b>	5	1X
<b>DI Water</b>	41	
<b>Tot Vol.</b>	50	

### After dPage:

<b>DNA sol.</b>	50	uL
<b>UV dilution factor</b>	30	
<b>Absorbance measured</b>	0.05535	
<b>Exc. Coeff.</b>	408303	
<b>Circular Strand</b>	4.06683272	uM

### Staple Strand Binding:

	<i>Volume (uL)</i>	<i>Conc. (uM)</i>
<b>Staple Strand</b>	15.78947368	24
<b>Tot Vol.</b>	65.78947368	

### Staple Strand Binding:

	<i>Volume (uL)</i>	<i>Conc. (uM)</i>
<b>Circular Strand</b>	10	0.81336654
<b>Staple Strand</b>	0.488019926	0.97603985
<b>Phi29 Buffer</b>	5	1x
<b>DI Water</b>	34.51198007	
<b>Tot Vol.</b>	50	

### Ligation:

	<i>Volume (uL)</i>	<i>Units</i>
<b>Ligase</b>	3	1200
<b>Tot Vol.</b>	68.7894737	



### Polymerisation:

	<i>Volume (uL)</i>	<i>Units</i>
<b>dNTPs</b>	2	
<b>Phi29</b>	1.5	15

# Characterization of RCA product

## Characterization:

- The designed template complementary yields G-Quadruplex for every one cycle of amplification.
- ThioflavinT (ThT) a fluorescent organic molecule is used for the characterization of G-Quadruplex.
- ThT has an excitation of 425 nm and an emission at 495 nm.
- In the presence of G-Quad the fluorescence of ThT gets enhanced by a large order and hence can be used as a quantifying agent.






**Spectrophotometer Results:**

Excitation: 425 nm  
Emission: 490 nm

Rough concentration estimate: 14 nM of DNA in cuvette  
Actual concentration: 308 nM  
Total Volume after RCA: 50  $\mu$ L  
Amount of DNA produced: 0.22  $\mu$ g  
\*Note: This calculation is only for G-Quad sequence



## Fluorescence measurements on G-quadruplex DNA



Excitation: 425 nm

Emission: 490 nm


Concentration of ThioflavinT in all samples: 0.5  $\mu$ M

Concentration's of DNA used: 0.5 nM, 1 nM, 2 nM, 5 nM, 8 nM, 10 nM, 15 nM and 20 nM

Buffer Used: TAEM 6mM Mg



# RCA reaction kinetics




## Details:

- RCA Initiation done in bulk and transferred to cuvette
- DNA amplification → Increase in G-quad → ThT binding to G-quad → Fl. increase

Excitation: 425 nm  
Emission: 495 nm  
Time monitored: 5 hrs  
Instrument: Spectrofluorometer


# RCA Summary:

Phosphorylation of Linear Strand  
Staple strand pairing  
Ligating the circular strand



Initial Conc.: 20  $\mu\text{M}$

Purification of circular strand (dpage)  
Ethanol concentration of DNA




**Critical Step:**


- Identifying band in gel.
- Extraction using electrodialysis.

Conc. after  
extraction:  $\sim 4.5 \mu\text{M}$

Staple strand pairing  
RCA



# RCA in GUV



## GUV Synthesis:

- Pipetting and centrifuging to make bilayer
- Time after reaction initiation = 10 min
- Flow in microchannel = 20 min

