

Last Update: 22 March 2023

2D Heterostructure Preparation




Jae-hoon Ji



Stamp Preparation

A. Prepare PDMS stamp




- Material (\$ 120)
 - PDMS material - Sylgard 184 (Dow Corning, USA), \$ 45
 - PPC material - Propylene carbonate (Sigma-Aldrich, CAS 25511-85-7) \$ 75
- Equipment (\$ 306, used for releasing bubbles inside PDMS to reduce surface roughness)
 - Vacuum Desiccator (Southern Labwaer) \$ 56
 - Oil free vacuum pump (Rocker) \$ 250



Stamp Preparation

A. Prepare PDMS stamp



- PDMS preparation


- Prepare the pre-polymer base without bubble (Under vacuum condition for 0.5 hour)
- Mix pre-polymer base and cross linking curing agent (10:1~15:1)
- Pour the mixed solution in the container
- Bubble releasing (Under vacuum condition for 1 hour)
- Cure the solution for 4 hour, 4 C
- If the surface of the stamp is not clean enough => Plasma treating

- PPC preparation


- Coat the surface of the stamp

Stamp Preparation

B. Attach the stamp on the glass





C. Transfer TMDC flakes on the stamp



Transfer method using micromanipulator

D. Transfer TMDC flakes on the substrate using micromanipulator and optical camera




Transfer Equipment

- Optical scope (AMScope, \$850)



Micropositioner




- up to 225x magnification
- Have to find fairly large flakes
- Ask the provider whether we can get more high magnification range.

Transfer Equipment

(a) Using Long Working Distance Objective

(a)



Objective – too expensive

+1	Qty	Docs	Part Number - Price	Available / Ships
	<input type="text"/>		MY5X-802 \$695.00	✓ 3-5 Days
	<input type="text"/>		MY10X-803 \$873.00	✓ Today
	<input type="text"/>		MY20X-804 \$2,056.00	Lead Time
	<input type="text"/>		MY100X-806 \$3,467.00	✓ 3-5 Days

Micropositioner

